

**Megan's PCR troubleshooting guide**  
**\*\*\*Suggestions are not in any particular order\*\*\***

**Case of the missing (previously present) PCR product**

- 1) Use new reagents  
dNTPs, primers, Taq are prone to degradation
- 2) Decrease annealing temp. by 4-8C
- 3) Increase or decrease template concentration
- 4) Increase primer concentration
- 5) Try different thermal cycle machine

**Product is longer than expected**

- 1) Length increased by intron
- 2) Wrong product being amplified (e.g. plant instead of fungus)

**Multiple bands (<50 bp or primer dimers)**

- 1) Use less primer, DNA template or Taq  
Add more DNA if primer dimers
- 2) Decrease annealing time
- 3) Increase annealing temp.
- 4) Decrease extension time

- 5) Decrease extension temp.
- 6) Increase KCl concentration
- 7) Increase MgCl<sub>2</sub> concentration
- 8) Add up to 5% DMSO

### **Multiple bands (>50 bp)**

- 1) Use less primer, DNA template or Taq
- 2) Increase annealing time
- 3) Increase annealing temp.
- 4) Increase extension time
- 5) Increase extension temp.
- 6) Decrease KCl concentration
- 7) Increase MgCl<sub>2</sub> concentration

### **Weak PCR product**

- 1) Optimize annealing temperature
- 2) Increase primer, DNA template, Taq concentrations (separately, to isolate problem)
- 3) Increase primer length (if oligo is very short)

### **PCR product smearing (on-line tips)**

- 1) Dilution series of template
- 2) Reduce number of cycles
- 3) Use new reagents and clean pipettes
- 4) Try suggestions for multiple bands >50 bp
- 5) Lower enzyme concentration
- 6) Optimize MgCL2 concentration

<b>Final Mg<sup>2+</sup> concentration in reaction (mM)</b>	<b>1</b>	<b>2</b>	<b>2</b>	<b>3</b>	<b>3</b>	<b>4</b>	<b>4</b>	<b>5</b>
	<b>5</b>	<b>0</b>	<b>5</b>	<b>0</b>	<b>5</b>	<b>0</b>	<b>5</b>	<b>0</b>
Required volume of 25 mM MgCl <sub>2</sub> per reaction (ul)	0	2	4	6	8	10	12	14

- 6) Optimize primer concentration or use fresh primers
- 7) Optimize primer design